



NovaPED sports Bike

Item no.
987686-000



With high-quality micro fibre cover



Forefoot padding

Advanced raised outer edge

Firm midfoot support

Heel stabilisation as a cup




NovaPED sports Cross-Country Skiing

Item no.
wide - 987657-000
narrow - 987656-000



With ARU-X-static cover, antibacterial



Flexible forefoot

Pressure reducing pad

Longitudinal arch support

Firm heel cup




NovaPED sports Golf

Item no.
wide - 987618-000
narrow - 987617-000



With high-quality micro fibre cover



Pressure reducing pad

Metatarsal support

Continuous padding layer

Arch support

Relief Basis V

Heel pad

Rear foot guide



NovaPED sports Downhill Skiing

Item no.
987635-000



With ARU-X-static cover, antibacterial



Forefoot padding

Metatarsal support

Longitudinal arch support

Relief Basis V

Heel padding




NovaPED sports Skating

Item no.
wide - 987632-000
narrow - 987631-000



With ARU-X-static cover, antibacterial



Continuous padding

Outer edge support

Longitudinal arch stabilisation

Raised cup




NovaPED sports Snowboard

Item no.
987638-000



With ARU-X-static cover, antibacterial



Pressure reducing pad

Metatarsal support

Arch supports

Very high cup



Schein Orthopädie Service KG
Hildegardstr. 5
42897 Remscheid, Germany
Tel. +49 2191 910-0
Fax +49 2191 910-100
info@schein.de
www.novapedsports.de



Germany since 1879

NovaPED sports
by schein



**Foot support technology
in your sports shoe**

Handmade in Germany

Specialist vendors

04.08.40.05.12.13
099993-201


Foot support technology in the sports shoe

The wide variety of sports and the different strains on the musculoskeletal system show how different the strain on our feet can be as well.

To avoid damage from excessive strain or alleviate existing complaints, we have developed the NovaPED sports foot support further. These are high-quality plastic blanks characterised by a low height.

The NovaPED sports foot support blanks are designed for different sports with their specific movement and strain patterns. They ensure high wearing comfort, prevent excessive strain and contribute to the alleviation of complaints.

How forces may act on the feet during sports and what movement patterns result is shown by these examples:



Movement analysis on the treadmill



Choosing the right sports shoe is important for the shoe and foot support to form a perfectly working unit. For a good fit in the shoe, the present insole must be taken out before the NovaPED sports foot support is adjusted to the shoe.

Analysing movement

A foot scan, an examination, a movement analysis and a detailed interview form the basis for individual foot support supply aligned with the customer's needs. They can be supplemented with a dynamic pressure measurement and video analysis.

The NovaPED sports foot support can be reworked for individual adjustment.

NovaPED sports Running

Item no. wide - 987620-000
narrow - 987624-000

With high-quality micro fibre cover

NovaPED sports Running with raised outer edge

Item no. wide - 987627-000
narrow - 987626-000

With high-quality micro fibre cover

NovaPED sports Running with cut-out of the plantar fascia tendon

Item no. 987623-000
Construction kit Item no. 987621-000

With high-quality micro fibre cover

NovaPED sports Football

Item no. 987610-000

With ARU-X-static cover, antibacterial

NovaPED sports Ball Sports

Item no. wide - 987661-000
narrow - 987668-000

With high-quality micro fibre cover

NovaPED sports Racket

Item no. breit - 987616-000
schmal - 987615-000

With high-quality micro fibre cover